


**REF: ME13-035**

English .....	1
Enska .....	6
EspaÑol .....	11
Dansk .....	16
Deutsch .....	21
Français .....	26
Italiano .....	31
Język polski .....	36
Nederlands .....	41
Norsk .....	46
Português .....	51
Suomi .....	56
Svenska .....	61

## ASSEMBLY INSTRUCTIONS VACUUMPUMP

## Instructions for use

### Low-pressure hose assembly for use with air for driving surgical tools

Important information on medical devices of Medex Instrumenten Service

## EXPLANATION OF SYMBOLS

These symbols will be used in the IFU, on the label of the device and/or on the device itself.

	Consult instructions for use		CE marking of conformity for medical devices
	Global Trade Item Number (GTIN) used by industry and health care for unique product identification. Also UDI-DI.		Distributor
	Catalogue number		Manufacturer
	Batch code		Date of manufacture
	Serial number		Non-sterile
	Unique Device Identifier		Not for general waste
	Medical Device		Use-by date
	Swiss Representative		Importer
	<p><b>Warning:</b> Indicates a potentially hazardous situation, which if not avoided, could result in death or serious injury. <b>Caution:</b> Indicates a potentially hazardous situation, which if not avoided, could result in minor or moderate injury. <b>Attention:</b> Indicates a potentially hazardous situation, which if not avoided, could result in property damage.</p>		

## GENERAL

**Non-observance of this document can lead to prolongation of the operation time. Read and understand all the safety and operating instructions contained in this IFU. If you do not understand these instructions, or have any questions contact your distributor before attempting to use this device.**

Improper handling and care as well as improper use can lead to premature wear of the medical device and/or to risks for patient and user.

It is the responsibility of the operator that all persons handling the device have understood and followed the instructions and directions in this document.

Be sure to read, understand and follow the instructions given below!

- In particular, be sure to heed all cautions, warnings and danger notices.
- Keep this document accessible to users at all times.
- This text applies equally to male and female persons. The double spelling has only been omitted for the sake of better readability.
- In this document, instrument(s) and instrumentation have the same meaning as medical device or medical devices.

## RECLAMATION AND WARRANTY

We thank you for having decided to buy a medical device from Medex Instrumenten Service, hereafter also referred to as MEDEX. This product carries the CE mark, which means that it satisfies the essential requirements laid down in the regulation (EU) 2017/745 for medical devices.

Medex Instrumenten Service is the manufacturer of this medical device.

**Medex Instrumenten Service**

Markerkant 10-136<sup>E</sup>  
1316 AL Almere  
Nederland  
Tel.nr. 0031 (0)36 534 12 98  
postbus@medex-instrumenten.com

The goods must be checked for completeness and potential damage in transit immediately upon receipt. Notice of any such damage must be given instantly.

Should you have any questions on how to handle the product or use it for clinical applications, please do not hesitate to contact your local distributor or MEDEX directly by phone or email:

Tel.: 0031 (0)36 534 12 98

Email: postbus@medex-instrumenten.com

The team at Medex Instrumenten Service does its utmost to ensure the delivery of high-quality medical devices. Nevertheless, if you have any comments and/or a complaint about our products, please do not hesitate to contact us using the above contact options, stating the article number (REF) and the lot number (LOT) and/or serial number (SN).

Any serious incident caused by the device must be reported to MEDEX and the competent authority of the Member State where the user and/or patient is/are located. Clearly state the classification of the incident and, where possible, the article number (REF) and the lot number (LOT) and/or serial number (SN) of the device.

The Terms and Conditions of the Metaalunie effective at the time shall apply. Agreements deviating from these Terms and Conditions do not restrict the legal rights of the buyer.

Any warranty exceeding the above provisions shall require a contractual form and shall exclude component-related vandalism and consumables.

MEDEX as manufacturer or distributor will not and cannot be held responsible for any direct, indirect, unforeseen and/or consequential damages due to injudicious and/or incompetent use and careless handling. MEDEX as manufacturer or distributor will not and cannot be held responsible for any direct, indirect, unforeseen or consequential damages due to incorrect handling, preparation, cleaning and maintenance of the medical devices.

If MEDEX devices are repaired by firms or persons without MEDEX authorisation for the repair, the warranty will become null and void. In the event of a change/modification of a MEDEX device, the CE and the warranty will become void and Medex Instrumenten Service cannot be held liable for any consequences/incidents that occur after change/modification of the MEDEX device.

Failure to comply with the instructions, improper handling or incorrect and inappropriate use of medical devices supplied and/or manufactured by MEDEX will result in exclusion of all warranty and/or liability claims.

## RECEIVING INSPECTION

The packaging of the device must be intact upon receipt. Remove product from package and inspect for damage. Verify that the model received is in working order. Damaged or incorrect products may not be used and must be returned to your local distributor or MEDEX. After thorough inspection of the product, please note any and all caution and warning labels in this manual.

- It is very important to allow the product to remain in original packaging for 24 hours to acclimatize to room temperature before use.
- When storing the hose assembly, ensure it is kept in a sealed package to avoid environmental damage. Store in a clean and dry environment to prevent damage.

## INTENDED USE

### Intended use

The low-pressure hose assemblies manufactured by MEDEX are intended for conveying (low-pressure) medical air, also known as instrument air or air for driving surgical tools, from wall socket to Zimmer Biomet Vacuum Pump REF 4232 and REF 4235 in operating rooms.

### Users

This device is to be used under the direction of a healthcare professional. The device may only be used effectively in the medical fields by adequately trained and qualified personnel. The attending physician or the user is responsible for the selection of the instruments for certain applications and the (operative) deployment, the adequate training and information and for the sufficient experience in handling the device.

### Environmental conditions

Use can take place under operating conditions in a conditioned OR according to national laws and regulations. Store the hose assembly in a clean and dry area.

These hose assemblies have not been evaluated for safety, heating, migration or compatibility in the magnetic resonance (MR) environment.

## USE

### Product description

The low-pressure hose assemblies are available for different wall sockets, but always in combination with a 5 meter hose and pneumatic fitting on the other side. For more information about the various models, please contact your local distributor or MEDEX.

### Product specifications

The low-pressure hose assemblies manufactured by MEDEX comply with EN-ISO 5359 and are made with a 5-layer medical grade PVC tubing containing an inner layer of antistatic PVC and reinforcing layer with polyester. The maximum working pressure is 200 psi ( $\pm 14$  Bar) at 21 °C (70 °F). The fittings are permanently crimped on each end with ferrules.

### Combination products

 Warning These hose assemblies are specifically manufactured to be used with the Zimmer Biomet Vacuum Pump REF 4232 and REF 4235. Always be sure to read, understand and follow the instructions given in the IFU of the Zimmer Biomet Vacuum Pump **before** the use of this device.

### Operating instructions

The attending physician or the user is responsible for the correct application of the device.

1. Inspect the hose assembly for damage or deviations from standard conditions. If any is found, do not use the hose assembly.
2. Connect the hose assembly to the correct mating wall socket for low-pressure medical air.
3. Connect the hose assembly to the vacuum pump until the second click. See page 'Assembly instructions Vacuumpump'.
4. Ensure that all connections are well assembled and free of leaks prior to use.

- Make sure the hose (and connector) do not get pressed between wall and inlet to avoid occlusion and/or deformation.
- Turn the NIST nut only by hand. DO NOT use a tool.

### Completion

After completion of the treatment/procedure:

1. To disconnect the hose assembly from the vacuum pump, first press the button to release excess air, press a second time to release the hose assembly from the vacuum pump.
2. Disconnect the hose assembly from the mating wall socket.
3. Put away the device in accordance with local regulations and/or hospital guidelines. Device should be cleaned after every use according to the cleaning instructions following in this document.

## **WARNING**

- Never use a device of which components are damaged or deviate from standard conditions. Breakage, cracks, damage, misuse or incorrect handling of components can cause the device to malfunction.
- Improper maintenance, handling, misuse and/or poor cleaning procedures can lead to premature wear and tear and can cause the device to become unsuitable for its intended use.
- This device may only be used by authorised skilled and qualified personnel according to specific use in the medical field.
- The user must ensure that the device functions as intended and is safe to use.
- Never use grease, oil, organic lubricants or flammable materials on or near the hose assembly.
- Never use any type of flame or flammable or explosive material near the hose assembly.
- Always follow EN-ISO 5359 for hose assembly and medical gas products.

## **PRECAUTIONARY MEASURES**

- Never use this device for activities for which it is not intended.
- This device should be closely examined for functionality or damage prior to use. A damaged device must not be used. Make sure a back-up device is available.
- Care must be taken to ensure that this instrumentation continues to function properly. The device must not be kinked or damaged in any way. Improper use or care of the device can cause the device to malfunction or not to function at all.
- This device should be cleaned and stored in a clean and dry environment after each use.

## **CLEANING INSTRUCTIONS**

Process stage	Relevant Aspect	Instruction
Before Processing	 Attention	<b>DO NOT:</b> <ul style="list-style-type: none"> <li>• Submerge product in any form of liquid. This will cause damage and void the warranty.</li> <li>• Use alcohol containing solutions.</li> <li>• Clean with pungent hydrocarbons.</li> <li>• Sterilise.</li> </ul>
Preparation before processing	Personal protective equipment	Use personal protective equipment (gloves, water repellent protective skirt)
Manual Cleaning and disinfection	Accessories	Lint-free wipes, soft bristle brushes, clean tips
	Cleaning and disinfectant agent	Either use lint-free wipes, prepare disinfectant solution (with cleaning effects) following the recommendations of the manufacturer and soak the wipes in the solution or use ready to use wipes already moistened with the disinfectant.
	Procedure	<ul style="list-style-type: none"> <li>• Thoroughly wipe the surfaces of the medical device with wipes sufficiently moistened with the disinfectant.</li> <li>• Always work from top to bottom and from clean to dirty areas.</li> <li>• Make sure that all surfaces, joints, holes, crevices etc. of the medical device can be reached during the cleaning/disinfection process and are completely wetted.</li> <li>• The surfaces shall stay wet for the time recommended for sufficient disinfection.</li> <li>• Treat areas which may not be reached with the wipe with a soft bristle brush and/or clean tips e.g. crevices and holes</li> <li>• Check for cleanliness, if soiling is still visible, repeat the above steps with new wipes.</li> <li>• Moisten sufficient lint free wipes with tap water to wipe surfaces for removal of disinfectant residues.</li> <li>• Let the medical device dry or use lint-free wipes to dry surfaces.</li> </ul>
Inspection	Visual inspection	Inspect visually the cleanliness for possible disinfectant residues

## INSPECTION AND MAINTENANCE

### Preventive maintenance

Devices should always be inspected by the user prior to treatment. Examination should be thorough and must include a visual and functional inspection. No leakage must occur. Never use any device with obvious signs of excessive wear or damage, or any device that is incomplete or otherwise not functioning.

### Yearly inspection

Inspect product once a year for damage or visual imperfections. Inspections must be performed by certified technicians.

### Troubleshooting

Contact your local distributor for technical support.

### Expiry date

It is recommended to replace the hose assembly within 5 (five) years from the date of which it was manufactured. Replacement is the sole responsibility of the user.

## PRODUCT DISPOSAL

### Packaging

The packaging can be disposed of via the old paper and household waste and/or plastic depending on hospital guidelines.

### Product disposal

*End of lifecycle:* Maintenance and circularity have been taken into account when designing this device for a sustainable life cycle. This makes it possible to recycle the (materials used in this) device to a large extent. At the end of its life, please return the device to your local distributor or MEDEX.

CAUTION: Instrumentation must be decontaminated prior to disposal.

### National regulations

National regulations and waste disposal guidelines must be observed for all waste disposal measures.

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