





David Mosimann founded Bien-Air in 1959. His son, Vincent, currently presides over the Bien-Air Medical Technologies Group.

BIEN-AIR

CREATOR OF EXCELLENCE

SINCE 1959

First electric motors
1975



ORL-E-92
The legendary pedal developed with Professor Fisch

1992



SWISS
MADE  QUALITY



3

1960

1970

1980

2010



2007

OSSEODUO
The high-speed multifunction control unit



2014

PM2
Bien-Air reinvent the handpiece. Ergonomic design and state-of-the-art technology

1981

RHEOSTAT
First milling system for ear surgery



1959

The Bien-Air adventure begins in Neuchâtel

Offering reliable and easy to use instruments to simplify the work of practitioners and constantly improve patient comfort has been Bien-Air's mission since its creation in 1959.

With over 60 years' experience in ENT surgery and dentistry, Bien-Air has won global recognition for the performance and reliability of its products.

Ergonomics, precision and durability are at the core of the development of every new product. Constantly in tune with the day-to-day requirements of practitioners, Bien-Air has moved beyond its previous innovation boundaries, constantly setting the bar higher.

A true culture of excellence, which fits perfectly with the tradition of Swiss-Made products from the renowned Watch Valley. Bien-Air supplies high quality products, which can only be guaranteed by using state-of-the-art technologies and manufacturing processes.

Bien-Air is a human-scale family business. Its close relationships with its customers and partners are one of its most valuable assets. This is why every individual receives personalised attention.

Our promise is to continue our innovative approach in order to provide you with advanced instruments and solutions. Making your work easier, every day.

PM2 HANDPIECES EXCLUSIVE DESIGN



4



SHAPING THE FUTURE OF HANDPIECES

Instruments to match your expertise!
The PM2 handpieces give you the perfect combination of new design and state-of-the-art technology. For your comfort and that of your patients.

PM2 line



Conventional instrument



Up to 51% lighter than previous generation

Internal irrigation, reduced size, enhanced cooling

Multiple adjustable lengths

360° directional irrigation

Perfect handling and balance


PM2: REDEFINING ERGONOMICS INTERNAL IRRIGATION

These handpieces are up to 51% lighter and 37% shorter than the previous generation. Handling your instrument requires less effort, leading to greater accuracy in your movements and less fatigue.

The Bien-Air PM2 connection offers a noticeable reduction in size and a great flexibility. Each PM2 handpiece is compatible with Bien-Air NANO, RAPIDO and BASCH-1 micromotors. The PM2 handpieces offer you high-level performance and reliability.

Our heavy duty handpieces feature a unique irrigation system. The irrigation hose connection is at the base of the instrument and is internally channeled through to the distal end.


This design keeps your handpiece from overheating during long procedures. Your handpiece will remain comfortable and cool throughout lengthy or demanding surgical case.



The internal irrigation system cools the handpiece while ensuring perfect handling and visibility

5

PM2 HANDPIECES AHEAD OF THE FIELD



6 PM2 handpieces allow
3 bur extension lengths

ADJUSTABLE BUR

The PM2 burs are fitted with additional retentive locking grooves at the base. This means that the bur can be extended to 3 different lengths, giving you more access to delicate areas while providing better visibility without compromising bur stability.


SECURE LOCKING SYSTEM

The PM2 range of burs has been specially designed to be used with PM2 handpieces and micromotors at speeds of up to 80,000 rpm.

The triangular coupling prevents any loss of drive during heavy use. All torque is transferred for greater performance.



OSSEOSTAP UNIQUE PRECISION



ULTIMATE CONTROL

The Bien-Air **OSSEOSTAP** is particularly notable for its light weight, precision, ease of handling and reliability. The **OSSEOSTAP**'s nose has a unique structure, it tapers down to the distal end which provides rigidity for accurate drilling. The small dimension of the distal end and the curve of the nose allows good visibility of the drilling site. The **OSSEOSTAP** system has been designed for micro-shaping during delicate stapedotomy and ossiculoplasty operations on the middle ear.

The **OSSEOSTAP** foot pedal is a compact, standalone component that operates on standard AA batteries. This allows the surgeon to maintain precise control of the **OSSEOSTAP** handpiece without the need of a control unit or power outlet. Compact and ergonomic, the **OSSEOSTAP** foot pedal is portable, making it a valuable alternative to controlling the **OSSEOSTAP** with the **OSSEODUO** console.


The **OSSEOSTAP** handpiece integrates Bien-Air brushless, sensorless micromotor. This powerful micromotor rotates at speeds ranging between 2,000 rpm and 12,000 rpm, smoothly and without vibrations. Despite their reduced size, **OSSEOSTAP** burs have perfect concentricity, which further enhances precision and comfort.

The **OSSEOSTAP** handpiece, its motor, cable and tools have been designed to withstand multiple wash cycles and autoclave sterilisation processes.



S120

ERGONOMICS AND PRECISION



The S120 system is offered with a wide choice of single use and reusable shaver blades.

S120 HANDPIECE

The S120 handpiece is light weight, precise and easy to handle. With its ergonomics and integrated irrigation, which allows blades to be changed without disconnecting the irrigation, the S120 is marked out by its ease of handling and balance. It ensures high-precision operations, even in areas that are difficult to access.

You will be particularly impressed by its versatility and simplicity.

SIMPLY EFFICIENT

The S120 handpiece has been designed around a single concept: to make the instrument simple to use, so that the surgeon can forget the tool and concentrate on the operation itself.

With its integrated irrigation, it is not necessary to disconnect the irrigation line to change the blade.

A wide choice of shaver blades and burs is available for the S120 handpiece. Whether single use or reusable, they will meet all your needs.

OSSEODUO

SIMPLE AND PRATICAL



RAPIDO MICROMOTOR

Max speed: 80,000 rpm
 Sound level:
 Weight: 115 g



NANO MICROMOTOR

Max speed: 80,000 rpm
 Sound level:
 Weight: 80 g

The RAPIDO is a powerful, ventilated micromotor with excellent torque. It provides significant cutting performance even in challenging surgical situations. The directional self-ventilation provides low operating temperatures, consistent performance and comfort during long procedures.

RAPIDO and NANO micromotors are compatible with the range of Bien-Air PM2 handpieces.

The NANO micromotor is compact, lightweight and extremely quiet. It provides speed, torque, ergonomics, and less fatigue at the end of each surgical procedure.

With its convenient size and versatility, the **OSSEODUO** paves the way for an optimal surgical experience.

Developed by and for surgeons, the **OSSEODUO** combines simplicity of use with performance to provide a comprehensive portfolio of micromotors and instruments. One **OSSEODUO** feature allows you to rapidly select any connected micromotor from the console or the foot pedal.

The **OSSEODUO** is a multi-functional power console specifically designed for the following surgical specialties:

- ENT Surgery
- Neurosurgery
- Spine Surgery

Offering independence to the surgeon, the pedal controls the micromotor's speed, rotation direction and activates irrigation.

Ultra-precision control of motors and handpieces

High-speed shaping up to 80,000 rpm

Integrated irrigation with adjustable flow rate between 8 and 135 ml/min



Two micromotors can be fitted to the **OSSEODUO** simultaneously. It has been designed to ensure perfect command of the RAPIDO and NANO micromotors, giving you access to the PM2 range of compact handpieces.

The entire Drill range is also compatible, as are the S120 Shaver handpiece and the **OSSEOSTAP** stapedotomy and ossiculoplasty handpiece.





A RATIONAL, HIGH-PERFORMANCE SOLUTION

Along with robust and stable, the **OSSEODOC** is very simple to use. It has been exclusively designed for the following surgical specialties:

- Otology surgery
- Laterall skull-based surgery

Its reliability and performance have made it a benchmark tool.

There are 4 motor speed ranges, a choice of rotational direction, 2 irrigation options and 4 peristaltic pump delivery rates. All functions are also available on the pedal.

Indicator lamps and audible signals indicate the functions that have been activated.

OSSEODOC COMPACT AND FUNDAMENTAL



BASCH-1: GUARANTEED POWER

The BASCH-1 delivers a substantial torque with a max speed of 40'000 rpm.

The BASCH-1, based on proven Bien-Air technology, is compatible with the **OSSEODOC** and the previous generation consoles, and allows to upgrade your system to new PM2 line handpieces.



BASCH MICROMOTOR	
Max speed:	40'000 rpm
Weight:	115 g

OSSEOUNO FOR NOSE SPECIALISTS

COMFORT AND PRECISION

The Bien-Air **OSSEOUNO** Shaver system offers performance along with a refined reduced size.

The shaver can be controlled in continuous rotation mode or oscillating rotation mode with a choice of 5 speed ranges. The **OSSEOUNO** also offers highly precise control of the rotation speed up to 8,000 rpm.

The **OSSEOUNO** is fitted with an integrated peristaltic pump. The pump delivery is easily adjustable from 15 to 130 ml/min, by pressing a key.

- Reduced dimensions
- Integrated irrigation with variable flow rate between 15 - 130 ml/min in 3 stages
- A choice of 5 speed ranges



Lightweight, precision S120 handpiece

Precise, easy blade positioning during the operation



S120 HANDPIECE

The S120 handpiece is light weight, precise and easy to handle. With its ergonomics and integrated irrigation, which allows blades to be changed without disconnecting the irrigation, the S120 is marked out by its ease of handling and balance. It ensures high-precision operations, even in areas that are difficult to access.

You will be particularly impressed by its versatility and simplicity.

PRODUCTS OVERVIEW



Set **OSSEOSTAP** + Control Unit
1700331-001



Set **OSSEOUNO*** + S120
1700300-001



Set **OSSEODUO***
S120 + RAPIDO 1700521-001
S120 + NANO 1700522-001
RAPIDO 1700523-001
NANO 1700524-001



Set **OSSEODOC***
BASCH-1 + pump 1700525-001
BASCH + pump 1700168-001
BASCH 1700169-001



OSSEOSTAP + Cable
1600685-001



S120 Shaver + Cable
1700302-001
S120 Shaver
1600614-001
S120 Cable
1600629-001



80K Micromotor + cable
1700263-001
80K Micromotor
1600518-001
80K/RAPIDO/PM PERFO cable
1600519-001



NANO Micromotor + cable
1700501-001
NANO Micromotor
1600777-001
NANO cable
1600865-001



RAPIDO Micromotor + cable
1700503-001
RAPIDO Micromotor
1600895-001
80K / RAPIDO / PM PERFO cable
1600519-001



BASCH-1 Micromotor + cable
1700502-001
BASCH-1 Micromotor
1600884-001
BASCH Cable
1600342-001



BASCH Micromotor + cable
1700074-001
BASCH Micromotor
1600076-001
BASCH cable
1600342-001

16

17

INTRA STRAIGHT OR ANGLED ATTACHMENT



PMAM 1122 70mm
1600207-001



PMAM 1122 95mm
1600336-001



PMAM 1122 125mm
1600337-001

Ø 2,35 mm
L 44,5 mm
L 70 mm
L 95 mm
L 125 mm



PMRM 1122 70mm
1600206-001



PM 1124 46mm
1600612-001



PMAM 1221 70mm
1600295-001

PM2
L 70 mm
L 95 mm
L 125 mm



PM2-S70HD
1600831-001



PM2-S70
1600767-001



PM2-A70HD
1600830-001



PM2-A70
1600765-001



PM2-A95
1600768-001



PM2-A125
1600769-001

PM2 STRAIGHT OR ANGLED ATTACHMENT



Aquacare
1600617-006



Lubrifiuid
1600064-006



Irrigation lines
1100037-010



Clips for irrigation
on BASCH and BASCH-1
10 pces / box
1300064-010



Clips for irrigation
and suction lines on S120 /
80K / NANO / RAPIDO
10 pces / box
1302844-010



Sterilization tray for:
PM2 burs, handpieces and micromotors 1600948-001
Blades and S120 handpiece 1600828-001
Burs and OSSEOSTAP handpiece 1600669-001
INTRA burs, handpieces and micromotors 1600477-001

*Each set includes 1 control unit, 1 foot pedal, 1 irrigation stand, 1 power cable, 10 irrigation lines, and 10 irrigation clips.

CONTROL UNITS

MICROMOTORS

PM2 – 50K* BURS RANGE

(*maximum usable speed for the bur: 50,000 rpm)



STAINLESS STEEL		DIAMOND		COARSE DIAMOND	
Non sterile – 1 pce / box		Sterile – 5 pces / box		Non sterile – 1 pce / box	
∅ mm	Reusable	∅ mm	Single use	∅ mm	Reusable
0.6	1100290-001	1100416-005	1100247-001	2.3	1100304-001
0.8	1100291-001	1100417-005	1100248-001	2.7	1100305-001
1.0	1100292-001	1100418-005	1100249-001	3.1	1100306-001
1.4	1100293-001	1100419-005	1100250-001	3.5	1100307-001
1.8	1100294-001	1100420-005	1100251-001	4.0	1100433-001
2.3	1100295-001	1100421-005	1100252-001	4.5	1100434-001
2.7	1100296-001	1100422-005	1100253-001	5.0	1100435-001
3.1	1100297-001	1100423-005	1100254-001	6.0	1100436-001
3.5	1100298-001	1100424-005	1100255-001	7.0	1100437-001
4.0	1100299-001	1100425-005	1100256-001		
4.5	1100300-001	1100426-005	1100257-001		
5.0	1100301-001	1100427-005	1100258-001		
6.0	1100302-001	1100428-005	1100259-001		
7.0	1100303-001	1100429-005	1100260-001		

CARBIDE		CROSSCUT		CONICAL CARBIDE	
Non sterile – 1 pce / box		Non sterile – 1 pce / box		Non sterile – 1 pce / box	
∅ mm	Reusable	∅ mm	Reusable	∅ mm	Single use
0.5	1100261-001	1.4	1100438-001	6.0	1100275-001
0.8	1100262-001	2.3	1100439-001	7.0	1100276-001
1.0	1100263-001	3.1	1100440-001		
1.4	1100264-001	4.0	1100441-001		
1.8	1100265-001	5.0	1100442-001		
2.3	1100266-001	6.0	1100443-001		
2.7	1100267-001				
3.1	1100268-001				
3.5	1100269-001				
4.0	1100270-001				
4.5	1100271-001				
5.0	1100272-001				
6.0	1100273-001				
7.0	1100274-001				



DIAMOND	
Non sterile – 1 pce / box	
∅ mm	Reusable
0.6	1100389-001
0.8	1100390-001
1.0	1100391-001
1.4	1100469-001
1.8	1100470-001
2.3	1100471-001

PM2 – 80K* BURS RANGE

(*maximum usable speed for the bur: 80,000 rpm)



STAINLESS STEEL		DIAMOND		CARBIDE		FULL CARBIDE	
Sterile – 5 pces / box		Sterile – 1 pce / box		Sterile – 1 pce / box		Sterile – 1 pce / box	
∅ mm	Single use	∅ mm	Single use	∅ mm	Single use	∅ mm	Single use
0.6	1100351-005	0.6	1100365-001	0.5	1100392-001	0.5	1100475-001
0.8	1100352-005	0.8	1100366-001	0.8	1100393-001	1.0	1100462-001
1.0	1100353-005	1.0	1100367-001	1.0	1100394-001	1.5	1100463-001
1.4	1100354-005	1.4	1100368-001	1.4	1100395-001	2.0	1100464-001
1.8	1100355-005	1.8	1100369-001	1.8	1100396-001	2.3	1100465-001
2.3	1100356-005	2.3	1100370-001	2.3	1100397-001	3.0	1100466-001
2.7	1100357-005	2.7	1100371-001	2.7	1100398-001	3.5	1100483-001
3.1	1100358-005	3.1	1100372-001	3.1	1100399-001	4.0	1100467-001
3.5	1100359-005	3.5	1100373-001	3.5	1100400-001	4.5	1100484-001
4.0	1100360-005	4.0	1100374-001	4.0	1100401-001	5.0	1100476-001
4.5	1100361-005	4.5	1100375-001	4.5	1100402-001	6.0	1100468-001
5.0	1100362-005	5.0	1100376-001	5.0	1100403-001	7.0	1100477-001
6.0	1100363-005	6.0	1100377-001	6.0	1100404-001		
7.0	1100364-005	7.0	1100378-001	7.0	1100405-001		

COARSE DIAMOND		EXTRA COARSE DIAMOND	
Sterile – 1 pce / box		Sterile – 1 pce / box	
∅ mm	Single use	∅ mm	Single use
2.3	1100379-001	4.5	1100381-001
3.1	1100380-001	6.0	1100382-001



STAINLESS STEEL		DIAMOND		COARSE DIAMOND		EXTRA COARSE DIAMOND	
Sterile – 1 pce / box							
∅ mm	Single use						
1.0	1100446-001	1.0	1100320-001	3.1	1100325-001	3.8	1100444-001
1.8	1100447-001	1.8	1100321-001	4.0	1100326-001	4.7	1100445-001
2.3	1100448-001	2.3	1100322-001	5.0	1100327-001	6.0	1100490-001
3.1	1100449-001	3.1	1100323-001	6.0	1100328-001		
4.0	1100450-001	4.0	1100324-001				
5.0	1100451-001						
6.0	1100452-001						



STAINLESS STEEL		DIAMOND		COARSE DIAMOND		EXTRA COARSE DIAMOND	
Sterile – 1 pce / box							
∅ mm	Single use						
1.8	1100453-001	1.0	1100308-001	3.1	1100313-001	3.8	1100316-001
2.3	1100454-001	1.8	1100309-001	4.0	1100314-001	4.7	1100317-001
3.1	1100455-001	2.3	1100310-001	5.0	1100315-001	6.0	1100491-001
4.0	1100456-001	3.1	1100311-001				
5.0	1100457-001	4.0	1100312-001				
6.0	1100458-001						

OSSEOSTAP BURS



OSSEOSTAP CARBIDE BURS

ø mm	1 pce / box sterile	Color coding
0.5	1501222-001	Violet
0.6	1501223-001	Blue
0.7	1501224-001	Green
0.8	1501225-001	Yellow
1.0	1501226-001	Orange
1.4	1501227-001	Grey
1.8	1501228-001	Brown
2.3	1501229-001	Red

Set comprising 6 sterile burs,
various diameters:
- 2 carbide burs
- 3 diamond burs
- 1 perforator
1700332-001



OSSEOSTAP DIAMOND BURS

ø mm	1 pce / box sterile	Color
0.6	1501230-001	Blue
0.7	1501231-001	Green
0.8	1501232-001	Yellow
1.0	1501233-001	Orange
1.2	1501428-001	Black
1.4	1501234-001	Grey
1.8	1501235-001	Brown
2.3	1501236-001	Red

OSSEOSTAP STAINLESS STEEL PERFORATOR

ø mm	1 pce / box sterile	Color
0.35	1501240-001	Black



Sterilization bur color coding table
5 pces / box
1304358-005

INTRA BURS



70 MM

STAINLESS STEEL

ø mm	Non sterile – 1 pce / box		Sterile – 5 pce / box	
	Reusable	Single-use	Reusable	Single-use
0.6	1100059-001	1100168-005		
0.8	1100060-001	1100169-005		
1.0	1100061-001	1100170-005		
1.4	1100062-001	1100171-005		
1.8	1100063-001	1100172-005		
2.3	1100064-001	1100173-005		
2.7	1100065-001	1100174-005		
3.1	1100066-001	1100175-005		
3.5	1100067-001	1100176-005		
4.0	1100068-001	1100177-005		
4.5	1100069-001	1100178-005		
5.0	1100070-001	1100179-005		
6.0	1100071-001	1100180-005		
7.0	1100072-001	1100181-005		

CARBIDE

ø mm	Non sterile – 1 pce / box		Sterile – 5 pce / box	
	Reusable	Single-use	Reusable	Single-use
0.5	1100014-001	1100152-005		
0.8	1100015-001	1100153-005		
1.0	1100016-001	1100154-005		
1.4	1100017-001	1100155-005		
1.8	1100018-001	1100156-005		
2.3	1100019-001	1100157-005		
2.7	1100020-001	1100158-005		
3.1	1100021-001	1100159-005		
3.5	1100022-001	1100160-005		
4.0	1100023-001	1100161-005		
4.5	1100024-001	1100162-005		
5.0	1100025-001	1100163-005		
6.0	1100026-001	1100164-005		
7.0	1100027-001	1100165-005		

DIAMOND

ø mm	Non sterile – 1 pce / box		Sterile – 5 pce / box	
	Reusable	Single-use	Reusable	Single-use
0.6	1100000-001	1100138-005		
0.8	1100001-001	1100139-005		
1.0	1100002-001	1100140-005		
1.4	1100003-001	1100141-005		
1.8	1100004-001	1100142-005		
2.3	1100005-001	1100143-005		
2.7	1100006-001	1100144-005		
3.1	1100007-001	1100145-005		
3.5	1100008-001	1100146-005		
4.0	1100009-001	1100147-005		
4.5	1100010-001	1100148-005		
5.0	1100011-001	1100149-005		
6.0	1100012-001	1100150-005		
7.0	1100013-001	1100151-005		

COARSE DIAMOND

ø mm	Non sterile – 1 pce / box
2.3	1100232-001
2.7	1100233-001
3.1	1100234-001
3.5	1100235-001
4.0	1100236-001
4.5	1100237-001
5.0	1100238-001
6.0	1100239-001
7.0	1100240-001

CROSSCUT

ø mm	Non sterile – 1 pce / box
1.4	1100226-001
2.3	1100227-001
3.1	1100228-001
4.0	1100229-001
5.0	1100230-001
6.0	1100231-001



44.5 MM

DIAMOND

Non sterile – 1 pce / box	
ø mm	Reusable
1.0	1100184-001

SHAVER BLADES



Double serrated Single serrated Smooth



Concave
Convex



SINGLE USE STRAIGHT SHAVER BLADES

	Double serrated	Single serrated	Smooth
Ø mm	5 pcses / box	5 pcses / box	5 pcses / box
3.0	1600538-005	1600541-005	1600544-005
3.5	1600539-005	1600542-005	1600545-005
4.0	1600540-005	1600543-005	1600546-005

Usable length 110 mm



REUSABLE STRAIGHT SHAVER BLADES

	Double serrated
Ø mm	1 pce / box
3.0	1600579-001
3.5	1600580-001
4.0	1600581-001

Usable length 100 mm – Can be dismantled for optimal cleaning



15° - SINGLE USE CURVED SHAVER BLADES - DOUBLE SERRATED

	Concave	Convex
Ø mm	5 pcses / box	5 pcses / box
4.0	1600570-005	1600571-005

Usable length 110 mm



40° - SINGLE USE CURVED SHAVER BLADES - DOUBLE SERRATED

	Concave	Convex
Ø mm	5 pcses / box	5 pcses / box
3.0	1600548-005	-
4.0	1600550-005	1600559-005

Usable length 110 mm



60° - SINGLE USE CURVED SHAVER BLADES - DOUBLE SERRATED

	Concave	Convex
Ø mm	5 pcses / box	5 pcses / box
3.0	1600573-005	-
4.0	1600574-005	1600575-005

Usable length 110 mm



90° - SINGLE USE CURVED SHAVER BLADES - DOUBLE SERRATED

	Concave	Convex
Ø mm	3 pcses / box	-
4.0	1600721-003	-

Usable length 120 mm



120° - SINGLE USE CURVED SHAVER BLADES - DOUBLE SERRATED

	Concave	Convex
Ø mm	3 pcses / box	-
4.0	1600722-003	-

Usable length 120 mm



ADENOID SHAVER BLADES



40° - SINGLE USE CURVED ADENOID SHAVER BLADE

Double serrated
Convex window orientation
External tube: Ø 4.0 mm
5 pcses / box
1600572-005

Usable length: 120 mm



LARYNGEAL SHAVER BLADES



SINGLE USE LARYNGEAL SHAVER BLADE

Concave window orientation
External tube: Ø 3.5 mm
1 pce / box
Usable length 220 mm: 1600568-001
Usable length 270 mm: 1600569-001

SHAVER BURS



15° - SINGLE USE CURVED STAINLESS STEEL DCR SHAVER BUR

Cutting bur: Ø 3.0 mm / External tube: Ø 3.5 mm
3 pcses / box
1600577-003

Usable length 110 mm



15° - SINGLE USE CURVED DIAMOND DCR SHAVER BUR

Cutting bur: Ø 3.0 mm / External tube: Ø 3.5 mm
3 pcses / box
1600650-003

Usable length 110 mm



SINGLE USE STRAIGHT CYLINDRICAL SHAVER BUR

Cutting bur: Ø 3.0 mm / External tube: Ø 3.5 mm
3 pcses / box
1600576-003

Usable length 120 mm



60° - SINGLE USE CYLINDRICAL SHAVER BUR

Cutting bur: Ø 3.0 mm / External tube: Ø 3.5 mm
3 pcses / box
1600578-003

Usable length 130 mm



S W I S S  M A D E